

**Baxter**

# Overbed Table - TA270

Designed to be  
simple yet sturdy  
and ergonomic.

- The **Autolift** system from Baxter may help reduce the risk of entrapment and overturning. With a simple push, the table top lifts and moves securely into place.
- The H and U-shaped base supports maximum stability and are compatible in all care environments.
  - H-shaped design features a reduced width front wheelbase.
  - U-shaped design is ideal for use with a variety of chairs.
- Locking dual-directional tilted table top creates a reading position.
- Plastic tray is designed for easy cleaning and water retention from glass spills.



High Pressure  
Laminate (HPL) Tray



Plastic Tray



U-shaped foot



H-shaped foot

## TECHNICAL SPECIFICATIONS

<b>Height adjustment</b>	HPL: 730 – 1080 mm (without rails) Plastic: 740 – 1090 mm (H-shaped foot)   752 – 1102 mm (U-shaped foot)
<b>Sizes</b>	Fixed models fixed table: 960 x 450 mm Tilting models tilting portion: 740 x 450 mm / fixed portion: 220 x 450 mm
<b>Weight</b>	18.5 kg
<b>Maximum load</b>	16 kg
<b>Material</b>	Surface area: HPL tray – High Pressure Laminate; Plastic tray – Polypropylene and UV resistance Undercarriage: Metallic grey epoxy coated steel (with plastic cover for H-shaped foot)
<b>Casters</b>	4 x 50 mm double-band with 2 brakes
<b>Tilting range</b>	Model TA270 TT max. 30° locked positions / 360° full rotational
<b>Colours</b>	See colours sheet for HPL tray / Grey (RAL 7035) for plastic tray

**For more information please contact your Baxter representative.**

This document is destined solely for use by healthcare professionals. Products shown in this brochure are intended for use with patients in departments of healthcare establishments. Baxter recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the products. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these products.

Furniture manufacturer: Hill-Rom S.A.S. – Z.I. du Talhouet, 56330 Pluvigner, France

Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products.

**[Hillrom.Baxter.eu](https://www.hillrom.com/baxter)**

Hill-Rom S.A.S.

Z.I. du Talhouet, 56330 Pluvigner, France

Baxter and Autolift are trademarks of Baxter International Inc. or its subsidiaries.

EMA-CS00-250096 (v1.0) 07/2025